




A Randomised Controlled Trial of Dilapan-S[®] versus Foley Balloon for Pre-induction Cervical Ripening (DILAFOL Trial)



Study design

A single center, randomised, open-label trial.
419 women were randomised from November 2016 to February 2018.

Principal Investigator

Antonio F. Saad, MD, University of Texas Medical Branch, Galveston, Texas.

Objective

Confirm DILAPAN-S non-inferiority to the Foley balloon for pre-induction cervical ripening at term.

Baseline characteristics


Nulliparous represented 41.9 % of subjects in DILAPAN-S group and 46.9 % in the Foley balloon group.

**DILAPAN-S is a highly effective cervical ripening agent,
with a vaginal delivery rate of 81%.**

**DILAPAN-S offered significantly higher maternal satisfaction over
Foley balloon during cervical ripening. Induced women were better able
to sleep, relax and perform daily activities without limitations.**

First round cervical ripening success rate

High rate of successful cervical ripening was comparable between the groups




In the DILAPAN-S group the outcome was achieved by leaving 5 dilators in place for 13 hours.

5
dilators


13
hours

The criteria for evaluation of successful ripening were primarily cervical dilation (>3cm) and effacement (>60%) rather than a complete Bishop score assessment.

Vaginal delivery rate



Caesarean section rate




The results indicated a trend towards more vaginal deliveries and a lower Caesarean section rate in DILAPAN-S group vs Folley balloon.

Vaginal delivery within 24 hours

Almost 50% of women delivered within 24 hours of DILAPAN-S insertion.*

*Data on file.



Safety

Low occurrence of adverse events confirmed DILAPAN-S favorable safety profile.

	DILAPAN-S	Foley Ballon
Vaginal bleeding*	3.1%	0.9%
Rupture of membranes*	0.9%	0.9%
Tachysystole*	0%	0%
Non-reassuring fetal status*	0.5%	1.4%
5 min Apgar < 7	0.5%	0.5%
Maternal infectious co-morbidity**	0%	0%


* During cervical ripening interval

** Related to device used

Maternal satisfaction

The graph demonstrates the percentage of women who were always or often able to achieve the desired activity during cervical ripening.

DILAPAN-S Foley balloon



DILAPAN-S offers greater maternal satisfaction during cervical ripening, compared with Foley balloon; induced women were significantly better able to sleep, relax and perform daily activities without limitations. Additionally, DILAPAN-S does not protrude from the introitus and does not need to be kept under tension.



DILAFOL trial - DILAPAN-S key facts

87 %

first round cervical
ripening success rate

81 %

vaginal
delivery rate

47 %

vaginal delivery
rate within 24 hours

0 %

serious maternal
and fetal complications



REFERENCE: Saad A. et al. A randomized controlled trial of DILAPAN-S vs Foley balloon for preinduction cervical ripening (DILAFOL trial), Am J Obstet Gynecol 220(3):275.e1-275.e9., March 2019

 **Dilapan-S**[®]

 **medicem**

DILAPAN-S is manufactured by company MEDICEM Technology s.r.o.,
Czech republic (www.medicem.com)



DS/EN/008/June22